

Intelligence • Efficiency




User-Friendly Control Panel

- Software shortcut function keys, allowing you to switch between patient information, observation and examination, examination and analysis interfaces with one click, comprehensively improving work efficiency and operating experience.
- High-performance "all-in-one" azimuth key offers sensitive control over magnification and focal length, enabling precise laser distance measurement with one click. This feature maintains lens stability during operation, ensuring consistent image quality.


Automatic Ranging Guidance

Using advanced laser ranging technology, it intelligently prompts the optimal working distance, allowing operations to be completed in a single step. This eliminates the need for repeated lens adjustments, enabling efficient image acquisition work and ensure high-quality results.



Intelligent Gesture Recognition

- The innovative gesture control system minimizes lens shake during image acquisition, ensuring precise imaging and quick focus. This technology reduce the risk of cross-infection by minimizing contact during use.




C9 UHD/C12 UHD
Video Colposcope

All-Embracing Perspective Revolutionary Clarity UHD Series

Video Colposcope




Global Headquarters:
Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen 518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com
U.S. and Canada inquiries:
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com
© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



UHD Series


Video Colposcope



Ultra HD imaging


4K Camera System Provides Insights into Lesion Details

C12 UHD is equipped with the latest SONY 4K CMOS technology, featuring over 8 million pixels and a resolution of 3840x2160. This ultra-high-definition field of view meticulously presents the details of tiny lesions and vascular morphology, capturing every moment of the changes in acetic acid test and vascular characteristics with exceptional clarity.



Ultra-Wide Color Gamut for True-to-Life Colors

BT.2020 wide color gamut display offers a broader range of color levels, resulting in richer and more accurate color observation. The ultra-high color reproduction ensures that colors are natural and realistic, faithfully restoring the colors of human tissues.



Medical-Grade 4K Display for an Enhanced Visual Experience



8MP high resolution for fine picture quality



Supports 5 color spaces and 4 color gamut selections, making color presentation more accurate




Supports 7 high-precision gamma correction modes, ensuring optimal image effects



Medical-grade monitor, meeting all safety and quality standards for medical devices

Multispectral Light Source Imaging

Multispectral light source imaging technology utilize a combination of white light and narrow-band light source to enhance the visibility of cervical epithelial blood vessels and tiny lesions. It leverages hemoglobin's strong absorption characteristics of blue light and green light to highlight microstructural changes in abnormal-shaped blood vessels and clearly define the boundaries of lesions.



By using white light along with blue and green narrow-band light, the multispectral filter imaging further highlights the contrast of blood vessels at different depths. This aids in identifying blood vessel differences in cervical intraepithelial lesions across various levels, providing a clearer understanding of lesion characteristics.

